


## 2006 Reaction of (R)-(-)-carvone with benzylamine in the presence of Montmorillonite K-10 to a Schiff's base



### Classification

#### Reaction types and substance classes

Reaction of the carbonyl group in ketones

ketone, amine, Schiff's base, natural product, acid catalyst

#### Work methods

removal of water by azeotropic distillation, stirring with magnetic stir bar, filtering, distilling under reduced pressure, evaporating with rotary evaporator, heating with oil bath

### Instruction (batch scale 100 mmol)

#### Equipment

250 mL round bottom flask, water separator, reflux condenser, heatable magnetic stirrer, magnetic stir bar, rotation evaporator, distillation apparatus, vacuum pump, oil bath

#### Substances

(R)-(-)-carvone (bp. 230 °C)	15.0 g (15.6 mL, 100 mmol)
benzylamine (bp. 185 °C)	11.8 g (12.0 mL, 110 mmol)
cyclohexane (bp. 81 °C)	170 mL
Montmorillonite K-10	3 g

#### Reaction

150 mL cyclohexane, 15.0 g (15.6 mL, 100 mmol) carvone, 11.8 g (12.0 mL, 110 mmol) benzylamine and 3.0 g montmorillonite K-10 are filled in a 250 mL round bottom flask with a water separator and reflux condenser. Under stirring with a magnetic stirrer the reaction mixture is heated to reflux until no more water is separated (3-4 hours).

#### Work up

After cooling down to room temperature the suspension is filtered through a folded filter paper into a 250 mL round bottom flask and the residue is washed with 20 mL cyclohexane. If

the filtrate is not clear, the filtration is repeated. The solvent is evaporated at a rotary evaporator. A yellow liquid remains as crude product.

Crude yield: 22.3 g; GC purity 90%

The crude product is transferred into a 50 mL round bottom flask and is fractional distilled under reduced pressure (0.1 hPa).

Yield: 16.8 g (70.2 mmol, 70%), light yellowish opalescent liquid; bp 128-130 °C (0.1 hPa), oil bath temperature up to 175 °C; GC purity 98% (see analytics), distillation residue: 2.20 g, yellow glutinous oil.

### **Waste management**

#### **Recycling**

The cyclohexane of the reaction solution is collected and redistilled.

Montmorillonite K-10 can be used again after drying.

#### **Waste Disposal**

<b>Waste</b>	<b>Disposal</b>
aqueous phase from water separator	solvent water mixtures, halogen free
distillation residue	organic solvents, halogen free
Montmorillonite K-10	solid waste, free from mercury

#### **Time**

5 hours

#### **Break**


After filtration of the montmorillonite K-10

#### **Degree of difficulty**

Medium

## Analytics

### Reaction monitoring using the amount of separated water



Start of water separation is defined as point 0 of the time scale.

The reaction time varies according to the speed of heating and distillation.

## GC

### GC-conditions

column: Zebtron ZB-1, length 15 m, internal diameter 0.25 mm, film 0.25  $\mu$ m,  
(Phenomenex, Torrance, CA, USA)

inlet: injector temperature 300 °C; split injection

carrier gas: He, pre-column pressure 100 kPa

oven: start temperature 50 °C (2 min), 8 °C/min to 200 °C (5 min), 8 °C/min to 250 °C (15 min).




detector: FID, 256 °C, H<sub>2</sub> 33.9 mL/min; synth. air 322 mL/min; make-up-gas N<sub>2</sub>



integration: integrator 4290 (Thermo Separation Products)

Percent concentration was calculated from peak areas.


Regarding the GCs on the following page:

Retention time (min)	Substance	Peak area %		
		crude product	fore-run	pure product (main fraction)
4.5	benzylamine	3.2	1.4	0
8.8	carvone	0.3	1.0	0
20.0	Schiff's base	90.0	92.9	97.8



**GC of the crude product****GC of the fore-run****GC of the pure product (mainfraction)**

**$^1\text{H}$  NMR spectrum of the crude product (250 MHz,  $\text{CDCl}_3$ )** **$^1\text{H}$  NMR spectrum of the pure product (250 MHz,  $\text{CDCl}_3$ )**

$\delta$ (ppm)	Multiplicity	Number of H	Assignment
1.77	s	3	$\text{CH}_3$
1.96	s	3	$\text{CH}_3$
2.0 - 2.6	m	4	$\text{CH}_2$ ring
2.90	m	1	tertiary H
4.69	m	2	$=\text{N}-\text{CH}_2-\text{Ph}$
4.80	m	2	$\text{CH}_2=\text{C}$
6.29	m	1	$-\text{CH}=\text{C}$ ring
7.2 - 7.5	m	5	$\text{CH}$ arene

**$^{13}\text{C}$  NMR spectrum of the pure product (250 MHz,  $\text{CDCl}_3$ )**

$\delta$ (ppm)	Assignment
18.4	$\text{CH}_3$
20.7	$\text{CH}_3$
30.7	$\text{CH}_2$ ring
32.1	$\text{CH}_2$ ring
41.7	$\text{CH} - \text{C}(\text{CH}_3) = \text{CH}_2$ ring
54.0	$\text{CH}_2 - \text{N} =$
110.0	$\text{CH}_2 = \text{C}$
126.4	$\text{CH}$ arene
127.5	$\text{CH}$ arene
128.3	$\text{CH}$ arene
134.2	$\text{CH} = \text{C}$ ring
135.8	$\text{CH} = \text{C}(\text{CH}_3) - \text{C}$ ring
140.9	$\text{C}_{\text{quart}}$ arene
148.1	$\text{CH}_2 = \text{C}(\text{CH}_3) - \text{C}$
167.5	$\text{C} = \text{N}$ ring
76.5-77.5	solvent

**IR spectrum of the crude product (film)****IR spectrum of the pure product (film)**

( $\text{cm}^{-1}$ )	Assignment
3080, 3070, 3030	C – H – valence, arene
2970, 2920	C – H – valence, alkene
1640, 1620	C = C – and C = N – valence
1580, 1500	C = C – valence, arene