


4025 Synthesis of 2-iodopropane from 2-propanol



Classification

Reaction types and substance classes

nucleophilic substitution
 iodoalkane, alcohol

Work methods

heating under reflux, stirring with magnetic stir bar, shaking out, extracting, filtering, draining of gases, working with wash bottles, distilling with fractionating column, heating with oil bath

Instruction (batch scale 10 mmol)

Equipment

10 mL round bottom flask, reflux condenser, 3 wash bottles, heatable magnetic stirrer, magnetic stir bar, separating funnel, distillation bridge, 10 cm Vigreux column, oil bath

Substances

2-propanol (dry) (bp 82 °C)	0.601 g (0.765 mL, 10.0 mmol)
phosphorus, red	0.124 g (4.00 mmol)
iodine	1.29 g (5.08 mmol)
aqueous NaOH solution (0.5 M)	
molecular sieve 3Å	
diluted aqueous NaHCO ₃ solution	
sodium sulfate for drying	

Reaction

The procedure must be performed in a hood. The apparatus consists of a 10 mL round bottom flask with a magnetic stir bar and a reflux condenser. The reflux condenser is connected to three wash bottles one after the other. The first one, directly connected to the reflux condenser, is filled with molecular sieve, the second one is empty, the third one contains a 0,05 M aqueous sodium hydroxide solution.

1.29 g (5.08 mmol) iodine, 0.124 g (4.00 mmol) red phosphorus and 0.601 g (0.765 mL, 10.0 mmol) dry 2-propanol is added to the reaction flask. With stirring the reaction mixture is heated to boiling with an oil bath and refluxed for two hours.

Work up

After the reaction the reflux condenser is replaced by a distillation bridge and the crude product is directly distilled from the reaction mixture at a boiling temperature of 80-90 °C. The distilled product is transferred into a separating funnel and washed with diluted aqueous NaHCO₃-solution. The organic phase is separated and dried over sodium sulfate. After filtration the crude product is distilled through a 10 cm Vigreux column.

Yield : 1.21 g (7.12 mmol, 71%); bp 88 °C, colourless liquid; $n_D^{20} = 1.495$

Comments

The use of molecular sieve as drying agent is recommended, otherwise the yield will be lower.

If there are small amounts of non reacted 2-propanol in the crude product it is no longer detectable after the washing procedure.

Waste management

Recycling

Molecular sieve can be regenerated.

Waste disposal

Waste	Disposal
aqueous phase from extraction	solvent water mixtures, containing halogen
residue from the first distillation	dissolve in water, neutralize with diluted NaOH, then: solvent water mixtures, containing halogen
residue from the second distillation	organic solvents, containing halogen
sodium sulfate	solid waste, free from mercury
contents of the third wash bottle	domestic waste water

Time

4 hours

Break

Before the first distillation and between the first and second distillation

Degree of difficulty

Easy

Instruction (batch scale 100 mmol)

Equipment

50 mL two-neck round bottom flask, Thielepape apparatus, reflux condenser, 3 wash bottles, heatable magnetic stirrer, magnetic stir bar, separating funnel, distillation bridge, 10 cm Vigreux column, oil bath

Substances

2-propanol (dry) (bp 82 °C)	6.01 g (7.65 mL, 100 mmol)
phosphorus, red	1.24 g (40.0 mmol)
iodine	12.9 g (50.8 mmol)
aqueous NaOH solution (0.5 M)	
molecular sieve 3Å	
diluted aqueous NaHCO ₃ solution	
sodium sulfate for drying	

Reaction

The reaction must be performed in a hood. The apparatus consists of a 50 mL two-neck round bottom flask with a magnetic stir bar and on top a Thielepape apparatus combined with a reflux condenser. The reflux condenser is connected to three wash bottles one after the other. The first one, directly connected to the reflux condenser, is filled with molecular sieve, the second one is empty, the third one contains a 0,05 M sodium hydroxide solution. The Thielepape apparatus has a bottom frit with porosity P-0 or an extraction cone from glass with a bottom frit with porosity P-0. As an alternative the Thielepape apparatus may have a refluxing device, which is covered with glass wool.

12.9 g (50.8 mmol) of iodine is transferred into the Thielepape apparatus. 1.24 g (40.0 mmol) red phosphorus and 6.01 g (7.65 mL, 100 mmol) dry 2-Propanol is added into the reaction flask and with stirring the mixture is heated to boiling with the oil bath. By refluxing of solvent iodine in the reflux condenser is dissolved and transferred continuously into the reaction mixture. The oil bath can be removed after a short time, because the release of heat during the exothermic reaction results in refluxing without additional heating. When the reaction is finished the mixture is refluxed for additional 30 minutes.

Work up

At this stage the reflux of the Thielepape apparatus is closed and the crude product is directly distilled into the Thielepape apparatus. The distilled product is taken from the apparatus and transferred into a separating funnel. It is washed with diluted aqueous NaHCO₃-solution. The organic phase is separated and dried over sodium sulfate. After filtration of the drying agent the crude product is distilled through a 10 cm Vigreux column.

Yield: 14.1 g (82.9 mmol, 83%.); bp 88°C, colourless liquid; $n_D^{20} = 1.495$

Comments

The use of molecular sieve as drying agent is recommended, otherwise yields are lower. The transference of solvent through the Thielepape apparatus should occur fast.

If there are small amounts of non reacted 2-propanol in the crude product it is no longer detectable after the washing procedure.

Waste management

Recycling

Molekular sieve can be regenerated.

Waste disposal

Waste	Disposal
aqueous phase from extraction	solvent water mixtures, containing halogen
residue from the first distillation	dissolve in water, neutralize with diluted NaOH, then: solvent water mixtures, containing halogen
residue from the second distillation	organic solvents, containing halogen
sodium sulfate	solid waste, free from mercury
contents of the third wash bottle	domestic waste water

Time

4 hours

Break

Before the first distillation and between the first and second distillation

Degree of difficulty

Easy

Instruction (batch scale 1 mol)

Equipment

500 mL two-neck round bottom flask, Thielepape apparatus, reflux condenser, 3 wash bottles, heatable magnetic stirrer, magnetic stir bar, separating funnel, distillation bridge, 10 cm Vigreux column, oil bath

Substances

2-propanol (dry) (bp 82 °C)	60.1 g (76.5 mL, 1.00 mol)
phosphorus , red	12.4 g (0.400 mol)
iodine	129 g (0.508 mol)
aqueous NaOH solution (0.5 M)	
molecular sieve 3Å	
diluted aqueous NaHCO ₃ solution	
sodium sulfate for drying	

Reaction

The procedure must be performed in a hood. The apparatus consists of a 500 mL two-neck round bottom flask and on top a Thielepape apparatus combined with a reflux condenser. The

reflux condenser is connected to three wash bottles one after the other. The first one, directly connected to the reflux condenser, is filled with molecular sieve, the second one is empty, the third one contains a 0.05 M sodium hydroxide solution. The Thielepape apparatus has a bottom frit with porosity P-0 or an extraction cone from glass with a bottom frit with porosity P-0. As an alternative the Thielepape apparatus may have a refluxing device, which is covered with glass wool.

29 g (0.508 mol) of iodine is added into the Thielepape-apparatus. 12.4 g (0.400 mol) red phosphorus and 60.1 g (76.5 mL, 1.00 mol) dry 2-propanol is added into the reaction flask. Then with stirring the mixture is heated to boiling with an oil bath. By the refluxing of the solvent iodine is dissolved and transferred continuously into the reaction flask. The oil bath can be removed after a short time, because the release of heat of the exothermic reaction results in refluxing without additional heating. When the reaction is finished the mixture is refluxed for additional 30 minutes.

Work up

At this stage the reflux of the Thielepape apparatus is closed and the crude product is directly distilled from the reaction mixture into the Thielepape apparatus. The distilled product is taken from the apparatus and transferred into a separating funnel. It is washed with diluted aqueous NaHCO₃ solution and dried over sodium sulfate. After filtration of the drying agent the crude product is distilled through a 10 cm Vigreux column.

Yield: 144 g (853 mmol, 85%); bp 88°C, colourless liquid; $n_D^{20} = 1.495$

Comments

The use of molecular sieve as drying agent is recommended, otherwise yields will be lower. It is important, that the transference of solvent through the Thielepape apparatus should occur fast.

If there are small amounts of non reacted 2-propanol in the crude product it is no longer detectable after the washing procedure.

Waste management

Recycling

Molekular sieve can be regenerated.

Waste disposal

Waste	Disposal
aqueous phase from extraction	solvent water mixtures, containing halogen
residue from the first distillation	dissolve in water, neutralize with diluted NaOH, then: solvent water mixtures, containing halogen
residue from the second distillation	organic solvents, containing halogen
sodium sulfate	solid waste, free from mercury
contents of the third wash bottle	domestic waste water

Time

5 hours

Break

Before the first distillation and between the first and second distillation

Degree of difficulty

Easy


Analytics**Reaction monitoring**

Stopping the reaction to take a sample lowers the final yield! Monitoring the reaction can be performed by IR spectroscopy.


Sample preparation:

1 mL of the reaction mixture is taken with a pipette and distilled with a micro distillation apparatus. The distilled product is used for IR spectroscopy without solvent.


The disappearance of the OH-band at $\approx 3400\text{ cm}^{-1}$ is an indication of the completeness of the reaction

 ^1H NMR spectrum of the pure product (500 MHz, CDCl_3)

δ (ppm)	Multiplicity	Number of H	Assignment
1.82	d	6	CH_3
4.26	m	1	I-CH

^{13}C NMR spectrum of the pure product (75.5 MHz, CDCl_3)

δ (ppm)	Assignment
31.13	CH_3
21.23	I-CH

IR spectrum of the pure product (film)

(cm^{-1})	Assignment
2983, 2917	C-H-valence, alkane
1455	C-H-deformation
1370	C-H-deformation, $-\text{CH}_3$ symm.